



Anti- Human AMH(aa 453-560)

Product Name	Anti- Human AMH (aa 453-560)
Catalog No	IT-010-011M2
Source	Mouse monoclonal antibody, clone 3D3F10
Isotype	IgG1
Formulation	Each vial contains 100µg purified IgG in 100µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing Human AMH protein (IT-010-010Ep).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with Human AMH (aa 453-560) as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

