

# Anti-Human Serum amyloid A-1 (SAA-1)



|                |  |
|----------------|--|
| Product Name   | Anti- Human Serum amyloid A-1 (SAA-1)  |
| Catalog No     | IT-010-012M7   |
| Source         | Mouse monoclonal antibody, clone 4E5   |
| Isotype        | IgG1   |
| Formulation    | Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.  |
| Immunogen      | DNA vaccine expressing Human Serum amyloid A-1 (SAA-1) (IT-010-012Ep).   |
| Applications   | ELISA, IP, IF etc.   |
| Epitope        | Conformational   |
| Specificity    | Reacts with Human Serum amyloid A-1 (SAA-1) as tested.   |
| Storage        | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures.   |

