



# Anti-NP(Ebola/Zaire)

Product Name	Anti-NP(Ebola/Zaire)
Catalog No	IT-014-006
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 0.2 mg purified IgG in 0.1 ml PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	DNA vaccine expressing the nucleoprotein of (Ebola/Zaire) (a.a.1-739) (Genebank No. AF086833).
Applications	ELISA(1:1000-1:100,000), IP, etc.
Specificity	Reacts with NP (Ebola/Zaire). Cross-react with NP from other subtypes.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

**ELISA Profile of IT-014-006**

