



Anti-GP(Ebola)-Pan

Product Name	Anti-GP(Ebola)-Pan
Catalog No	IT-014-007
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 0.2 mg purified IgG in 0.1 ml PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	DNA vaccine expressing glycoprotein of (Ebola/Zaire) (a.a. 1-676) (GB# U23187) and (Ebola/Cote d'Ivoire) (a.a. 1-672) (GB# NC_014372).
Applications	ELISA (1:1000-1:100,000), WB (1:200-1:2000), IP, etc.
Specificity	Reacts with GP(Ebola/Zaire) and GP(Ebola/Cote d'Ivoire) glycoprotein. Cross-reactivity with other subtypes not determined.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing.</u> The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

