



VP40 (Ebola/Zaire)

Product Name	VP40 (Ebola/Zaire)
Catalog No	IT-014-011Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of protein in PBS with 8M Urea.
Protein	C-terminal 6xHis tagged Matrix protein VP40 (Ebola/Zaire) (GenBank # Q77DJ6, aa 1-326)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>90% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

