

Glycoprotein (Rift Valley Fever Virus)



Product Name	Glycoprotein (Rift Valley Fever Virus)
Catalog No	IT-017-001Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of protein in PBS with 8M Urea.
Protein	C-terminal 6xHis tagged Glycoprotein (RVFV) (a.a. 154-583 of the M protein. This region is also called Phlebovirus_G1) (Genebank Accession No: AAA47449).
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>85% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

