

Glycoprotein G2 (SFTS Virus)



Product Name	Glycoprotein G2 (SFTS Virus)
Catalog No	IT-017-005p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 50 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged Phlebovirus glycoprotein G2 (SFTS virus HN6) (a.a.563-1035) (Genebank No. ADZ04490).
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Glycoprotein G2 (SFTS Virus)