

Anti-NP(SFTS Virus)(Huaiyangshan Virus)



Product Name	Anti-NP(SFTS Virus)(Huaiyangshan Virus)
Catalog No	IT-017-006M11
Source	Mouse monoclonal antibody, clone 40C9
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing NP (SFTS Virus)(Huaiyangshan Virus) protein (a.a. 1-245) (Genebank Accession No. AEA29891; or ADZ04472).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts specifically with NP (SFTS Virus)(Huaiyangshan Virus). Minimal cross-reactivity to other Bunyavirus nucleocapsid protein.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-017-006M11

