

# Glycoprotein (Rift Valley Fever Virus)



Product Name	Glycoprotein (Rift Valley Fever Virus)
Catalog No	IT-017-008p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged Glycoprotein (Rift Valley Fever Virus) protein (a.a. 691-1139), also called "Phlebovirus_G2" (Genebank Accession No. AAA47449)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>90%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Glycoprotein (Rift Valley Fever Virus)