


E1 glycoprotein (Mayaro virus)



Product Name	E1 glycoprotein (Mayaro virus)
Catalog No	IT-022-024Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein in PBS with 8M Urea. Reconstitute the protein with 100 µl of Millipore water.
Protein	C-terminal V5-6xHis tagged E1 glycoprotein (Mayaro virus) antigen (a.a.1-413) (Genbank Accession No. NP_740694).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: E1 glycoprotein (Mayaro virus)