

# Anti- Human HER3 / ErbB3



---

Product Name	Anti- Human HER3 / ErbB3
Catalog No	IT-300-005M1
Source	Mouse monoclonal antibody, clone 6E4
Isotype	IgG1
Formulation	Each vial contains 100µg purified IgG in 100µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing Human HER3 / ErbB3 Protein (IT-300-005p).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with Human HER3 / ErbB3 Protein as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

