

Anti- Human c-KIT / CD117



| | |
|----------------|--|
| Product Name | Anti- Human c-KIT / CD117 |
| Catalog No | IT-300-019M2 |
| Source | Mouse monoclonal antibody, clone 9A3 |
| Isotype | IgG2a |
| Formulation | Each vial contains 100µg purified IgG in 100µl PBS with 40% glycerol |
| Immunogen | DNA vaccine expressing Human c-KIT / CD117 Protein (IT-300-019p). |
| Applications | WB, ELISA, IP, IF etc. |
| Epitope | Linear |
| Specificity | Reacts with Human c-KIT / CD117 Protein as tested. |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |

