

# tpn39b (T. pallidum/Nichols strain)



---

Product Name	tpn39b (T. pallidum/Nichols strain)
Catalog No	IT-501-010Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl with 8M Urea.
Protein	6xHis tagged tpn39b (Treponema pallidum/Nichols strain) (a.a. 1-361) (GenBank Accession No. AAC65966).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

