

Anti-Nucleocapsid protein N(CDV)



Product Name	Anti-Nucleocapsid protein N(CDV)
Catalog No	IT-007-001M5
Source	Mouse monoclonal antibody, clone 5E5D10
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Nucleocapsid protein N(CDV) (IT-007-001Ep).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with Nucleocapsid protein N(CDV) as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

