



## Anti-MB (Ureaplasma parvum)

Product Name	Anti-MB (Ureaplasma urealyticum)
Catalog No	IT-503-016M3
Source	Mouse monoclonal antibody, clone 2G7D6
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing MB (Ureaplasma parvum) (IT-503-016Ep).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts strongly with both MB (Ureaplasma parvum) and MB (Ureaplasma urealyticum) as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

