

Anti- Gc(CCHFV)



Product Name	Anti-Gc(CCHFV)
Catalog No	IT-017-020M1
Source	Mouse monoclonal antibody, clone 4F3B3
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Gc(CCHFV) (IT-017-020Ep).
Applications	ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with Gc(CCHFV) protein as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

