




HA(Δ TM)(B/Washington/02/2019)

Product Name	HA(Δ TM)(B/Washington/02/2019)
Catalog No	IT-003-B22 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	HA(Δ TM)(B/Washington/02/2019) (aa 18-545; see seq below) with a 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	\geq 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM)(B/Washington/02/2019)

Note:

Protein Sequence:

MKAIIVLLMVVTSNADRICTGITSSNSPHVVKATATQGEVNVTGVIPLTTTPTKSHFANLKGTTETRGKLCPCCLN
CTDLVALGRPKCTGKIP SARVSILHEVRPVTSGCFPI MHDRTKIRQLPNLLRGYEHVRLSTHNVINAEDAPG
RPYEIGTSGSCP NITNGNGFFATMAWAVPKNKATNPLTIEVPYICTEGEDQITVWGFHSDnETQMAKLYGD
SKPQKFTSSANGVTTHYVSQIGGFPNQTEDGGLPQSGRIVVDYMVQKSGKTGTITYQRGILLPQKVVWCSG
RSKVIKGSPLIGEADCLHEKYGGLNKS KPYTGEHAKAIGNCPIWVKTP LKLANGTKYRPPAKLLKERGFFGA
IAGFLEGGWEGMIAGWHGYTSHGAHGVAADLKSTQEAINKITKNLNSLSELEVKNLQRLSGAMDELHNEI
LELDEKVDLDRADTISSQIELAVLLSNEGIINSEDEHLLALERK LKMLGPSAVEIGNGCFETKHKCNQTCLDRI
AAGTFDAGEFSLPTFDSL NITAASLND DGLDNHT