

## Human Interleukin-36 gamma / IL36G



Product Name	Human Interleukin-36 gamma / IL36G
Catalog No	IT-300-163Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg protein in PBS.
Protein	Human Interleukin-36 gamma / IL36G (a.a. 18-169) (GenBank Accession No. Q9NZH8). No tag.
Applications	WB standard, antibody ELISA, immunogen, Functional assay, etc.
Purity	>95% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human Interleukin-36 gamma / IL36G