

Anti-Spike RBD Domain of SARS-CoV-2 (CoVID-19/Wuhan), Mouse Monoclonal



Product Name	Anti-Spike RBD Domain of SARS-CoV-2 (CoVID-19/Wuhan)
Catalog No	IT-002-031M1
Source	Mouse Monoclonal Antibody
Isotype	Mouse IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Plasmid vector expressing Spike protein of SARS-CoV-2 (CoVID-19/Wuhan).
Applications	ELISA
Specificity	Reacts with Spike Protein RBD Domain from Coronavirus SARS-CoV-2 (CoVID-19/Wuhan) (IT-002-031p).
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

