

Anti- Nucleocapsid (CoVID-19/Wuhan)



Product Name	Anti- Nucleocapsid (CoVID-19/Wuhan)
Catalog No	IT-002-033M1
Source	Mouse monoclonal antibody, clone 6B5B4
Isotype	IgG2b
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Nucleocapsid (CoVID-19/Wuhan) (IT-003-033Ep).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with Nucleocapsid (CoVID-19/Wuhan) as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

