


HA1(A/Hong Kong/45/2019)(H3N2)



| | |
|----------------|--|
| Product Name | HA1(A/Hong Kong/45/2019)(H3N2) |
| Catalog No | IT-003-00439p |
| Source | Recombinant protein purified from 293 cell culture |
| Formulation | Each vial contains 100 µg of purified protein (1 mg/ml) in PBS. |
| Protein | HA1(A/Hong Kong/45/2019)(H3N2)(aa 17-346; see seq below) with a 6xHis tag at C-terminus. |
| Applications | WB standard, antibody ELISA, Functional analysis, etc |
| Purity | ≥95% |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |



AA SEQ:

QKIPGNDNSTATLCLGHHAVPNGTIVKTTITNDRIEVTNATELVQNSSIGEICDSPHQILDGGNCTLIDALLGDPQ
CDGFQNKKWDLFVERSRAYSNCYPYDVPDYASLRSLVASSGTLEFKNESFNWAGVTQNGKSFSCIRGSSSSFF
SRLNWLTHLNYTYPALNVTMPNKEQFDKLYIWGVHHPGTDKQISLYAQSSGRITVSTKRSQQA VIPNIGSRP
RIRDIPSRISYWTIVKPGDILLINSTGNLIAPRGYFKIRSGKSSIMRSDAPIGKCKSECITPNGSIPNDKPFQNVNRI
TYGACPRYVKQNTLKLATGMRNVPEKQTRG