


HA1(A/Darwin/6/2021)(H3N2)



Product Name	HA1(A/Darwin/6/2021)(H3N2)
Catalog No	IT-003-00443p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	HA1(A/Darwin/6/2021)(H3N2)(aa 17-345; see seq below) with a 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



AA SEQ:

QKIPGNDNSTATLCLGHHA V P N G T I V K T I T N D R I E V T N A T E L V Q N S S I G E I C G S P H Q I L D G G N C T L I D A L L G D P Q C D G F Q N K E W
D L F V E R S R A N S N C Y P Y D V P D Y A S L R S L V A S S G T L E F K N E S F N W T G V K Q N G T S S A C I R G S S S S F F S R L N W L T S L N N I Y P A Q N V T
M P N K E Q F D K L Y I W G V H H P D T D K N Q I S L F A Q S S G R I T V S T K R S Q Q A V I P N I G F R P R I R G I P S R I S I Y W T I V K P G D I L L I N S T G N L I A P
R G Y F K I R S G K S S I M R S D A P I G K C K S E C I T P N G S I P N D K P F Q N V N R I T Y G A C P R Y V K Q S T L K L A T G M R N V P E K Q T R