

Anti-MOMP (C. trachomatis strain C/TW-3 Serovar C)



Product Name	Anti-MOMP (C. trachomatis strain C/TW-3 Serovar C)
Catalog No	IT-502-002M1
Source	Mouse monoclonal antibody, clone 3C4
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing MOMP (C. trachomatis strain C/TW-3 Serovar C) (IT-502-002Ep).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with MOMP (C. trachomatis strain C/TW-3 Serovar C) protein only as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

