

IMV surface protein A29L [Monkeypox virus MPXV]



Product Name	IMV surface protein A29L [Monkeypox virus MPXV]
Catalog No	IT-012-019Ep
Source	Recombinant protein purified from E.coli.
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl.
Protein	IMV surface protein A29L [Monkeypox virus MPXV] (aa 1-110, Genbank No. QJQ40281) with a 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

