

Anti-H3(WR101)(Vaccinia Virus)



Product Name	Anti-H3(WR101)(Vaccinia Virus)
Catalog No	IT-012-005M2
Source	Mouse monoclonal antibody, clone 1B6F10
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant H3(WR101)(Vaccinia Virus) (IT-012-005Ep).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with H3(WR101)(Vaccinia Virus).
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

