



p15 (T. pallidum/Nichols strain)

Product Name	p15 (T. pallidum/Nichols strain)
Catalog No	IT-501-001Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0) and 150mM NaCl.
Protein	p15 (T. pallidum/Nichols strain) (a.a. 18-141; GenBank# P16055) with N-terminal GST and C-terminal 6xHis tag.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: p15 (T. pallidum/Nichols strain)

Note: