

Hemagglutinin Protein H (CDV)



Product Name	Hemagglutinin Protein H (CDV)
Catalog No	IT-007-004p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	N-terminal 6xHis tagged Human Fc-Hemagglutinin Protein H (Canine distemper virus/strain Onderstepoort)(a.a.144-604)(Gene Bank Accession No.AAK54669) fusion protein.
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95%, based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

