

Polyclonal Antibody Production Services



Features and Advantages

- DNA vaccine
- No protein purification
- High titer

Applications and Service Scale

IT offers customized polyclonal antibody manufacture services. Antibody is generated using DNA vaccine process through IT's unique and proprietary expression system. Polyclonal antibodies generated by IT are suitable for:

- Western Blotting
- ELISA
- Other research applications

Order Information

- **Cost**
 - \$2000
- **Animal**
 - Two rabbits (New Zealand)
- **Duration**
 - 3-4 months
- **Product**
 - 5 ml pre-sera/animal
 - 0.5 ml 1st bleed sera after 3rd vaccination
 - 25-30 ml terminal bleed sera after 5 vaccinations.
- **Optional (\$500 extra)**
 - IgG purification.

Guarantee: if the production is not successful or the titer of the anti-sera is unreasonably low, we will refund you 50% of the cost.